

# BOAR BETTER®

PATENTED  
TRIPLE  
PHEROMONE  
BSA\*\*

## Facilitates reproductive performance and efficiency, especially when used with a boar.


- Works through olfactory receptors to induce sows estrus behavioral responses.
- Supplies boar pheromones that can enhance boar performance when the environment may not have effective boars available.
- Can replace the need for a live boar.
- Has been used safely on thousands of sows under commercial conditions with no negative health outcomes.

### COMPOSITION:

Swine pheromones BSA\* (quinolone, androstenone and androstanol)... 0.006%  
Inert ingredients (contains IPA) ..... 99.994%

½ Liter (16.9 oz)

\*\*BSA (Boar Saliva Analog)



**BOARBETTER** spray is a management tool for heat detection in breeding female swine. This product contains a patented solution of three pheromones blended in a unique combination referred to as Boar Saliva Analogue (BSA). These three pheromones occur naturally in the Boar. Boar pheromones are responsible for the "boar odor" and this product can reduce the need to maintain live "teaser" boars.

**BOARBETTER** was developed to mimic the three pheromones found in boars (quinolone, androstenone and androstanol). This patented combination can maximize the breeding female's mating behavior response.

Laboratory and field studies confirm the combination of compounds in BSA results in a synergistic benefit. **BOARBETTER** results in a greater lordosis behavior (reproductive behavior) in sows compared to single pheromone formulations.

IN PARTNERSHIP WITH  
TEXAS TECH UNIVERSITY



### DIRECTIONS FOR USE:

- Begin using **BOARBETTER** post weaning (recommended on day 3-4 after weaning). Applying earlier post weaning may reduce the wean-to-estrus interval.
- Spray 4 mL per application of **BOARBETTER** onto the female's snout (directing towards the nostrils) prior to heat checking.
- Apply back pressure test to check for signs of standing estrus.
- If using conventional artificial insemination breed the female if she is in standing heat.
- If using PCAI (post cervical artificial insemination) follow farm specific PCAI breeding protocols.
- If female is not in estrus, repeat treatment until estrus occurs.

**WARNING:** Not for use in humans. Keep out of reach of children. Do not ingest. Do not spray in human eyes, mouth or face.

**HANDLING:** Wear protective gloves, protective clothing, eye protection and face protection. Wash hands after each use. Keep away from heat/sparks, open flames/hot surfaces. – No smoking.

**STORAGE:** Store at room temperature in a well-ventilated place. Do not freeze. Do not expose to heat above 110° F.

**DISPOSAL:** comply with environmental regulations. Contact a licensed waste disposal service to dispose of this product.

\*Boar Saliva Analog – patented US Patent #US9480689, US2018/065171 \* Trial #AB2015-2

**MANUFACTURED FOR:**  
Animal Biotech  
1601 Elm Street, Suite 3500,  
Dallas, TX 75201



Questions or comments?  
Please contact us at [boar-better.com](http://boar-better.com)