

VET ONE

Distributed by MWI Animal Health

3041 W. PASADENA DR., BOISE, ID, 83705

Customer Service: 888-694-8381

Toll-Free: 800-824-3703

Website: www.mwianimalhealth.com

Website: www.vetone.net



Every effort has been made to ensure the accuracy of the information published. However, it remains the responsibility of the readers to familiarize themselves with the product information contained on the USA product label or package insert.

STERILE SALINE



VetOne

Solution

Sterile - Preservative Free

CAUTION: Federal law restricts this drug to use by or on the order of a licensed veterinarian.

WARNING:

KEEP OUT OF REACH OF CHILDREN

FOR ANIMAL USE ONLY

COMPOSITION: Each 100 mL of sterile aqueous solution contains:

Sodium Chloride	0.9 g
-----------------	-------

Milliequivalents per liter

Cations

Sodium	154 mEq/L
--------	-----------

Anions

Chloride	154 mEq/L
----------	-----------

Total osmolarity is 308 milliosmoles per liter.

CAUTION: This product contains no preservatives. Use entire contents when first opened. Discard any unused solution.

INDICATIONS: For use in replacement therapy of sodium, chloride and water which may become depleted in many diseases. Because this solution is isotonic with body fluids, it may also be used as a solvent or diluent for antibiotics and other pharmaceuticals and biologicals where compatible, and for washing mucous membranes and other tissue surfaces.

DOSAGE AND ADMINISTRATION: Warm to body temperature and administer slowly by intravenous or subcutaneous injection. The amount and rate of administration must be judged by the veterinarian in relation to the condition being treated and the clinical response of the animal, being careful to avoid overhydration. When used as a solvent or diluent for pharmaceuticals and biologicals, follow manufacturer's directions.

Store between 15°C and 30°C (59°F-86°F).

Distributed by: MWI, Boise, ID 83705

(888) 694-8381

www.vetone.net

Manufactured by: Bimeda-MTC Animal Health Inc., Cambridge, Ontario, Canada N3C 2W4

Net Contents:	
1000 mL	Iss.01.07

	8XXXXXX V1 501032
--	----------------------

CPN: 1315019.1